




Recombinant Chicken Anosmin-1 (ANOS1), partial

Product Code	CSB-EP011978CH
Relevance	May be an adhesion-like molecule with anti-protease activity.
Abbreviation	Recombinant Chicken ANOS1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33005
Alias	Kallmann syndrome protein homolog
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SPAGPGAATARRQDEAFSTARCTSRCLSLQITRISAFFKHFQNNGLAWCQN HKQCSKCLEPCKESWDLKKNHCQSFCEPLFPKKNYECLTSCEFLKYILSVKQG DCPAPEKASGFAAACVESCEADSECSGVKKCCSNGCGHTCQVPKNLYKGVP LKPRKELKFIELQSGDLEVKWSSKFNISIEPVIYVVQRRWNQGIHPSEDDATNW QTVAQTDDERVQLSDIRASRWYQFRVAAVNVHGTGRGFTAPSKHFRSSKDP
Research Area	Others
Source	E.coli
Target Names	ANOS1
Protein Names	Recommended name: Anosmin-1 Alternative name(s): Kallmann syndrome protein homolog
Expression Region	22-281aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	33.2kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.