




Recombinant Human Metallothionein-3 (MT3)

Product Code	CSB-EP015122HU
Relevance	Binds heavy metals. Contains three zinc and three copper atoms per polypeptide chain and only a negligible amount of cadmium. Inhibits survival and neurite formation of cortical neurons in vitro.
Abbreviation	Recombinant Human MT3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25713
Alias	GIFB Short name: GIF Growth inhibitory factor Metallothionein-III Short name: MT-III
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MDPETCPGPSGGSGCTCADSCKCEGCKCTSCKKSCCSCCPAECEKCAKDCVC KGGEAAEAEAEKCSCCQ
Research Area	Neuroscience
Source	E.coli
Target Names	MT3
Protein Names	Recommended name: Metallothionein-3 Short name= MT-3 Alternative name(s): GIFB Short name= GIF Growth inhibitory factor Metallothionein-III Short name= MT-III
Expression Region	1-68aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	22.9kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.