




Recombinant Human Tropomodulin-1 (TMOD1)

Product Code	CSB-EP023909HU
Relevance	Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton. May play an important role in regulating the organization of actin filaments by preferentially binding to a specific tropomyosin isoform at its N-terminus.
Abbreviation	Recombinant Human TMOD1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28289
Alias	Erythrocyte tropomodulin
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSYRRELEKYRDLDEDEILGALTEEELRTLLENELDELDPDNALLPAGLRQKDQ TTKAPTGPFKREELLDHLEKQAKEFKDREDLVPYTGEKRGKVWVWPQKPLDP VLESVTLEPELEEALANASDAELCDIAAILGMHTLMSNQQYYQALSSSSIMNKE GLNSVIKPTQYKVPDDEPNSTDVEETLERIKNNDPKLEEVNLLNIRNIPIPTLKA YAEALKENSYVKKFSIVGTRSNDPVAYALAEMLKENVLKLTLNVESNFISGAGIL RLVEALPYNTSLVEMKIDNQSQPLGNKVEMEIVSMLEKNATLLKFGYHFTQQG PRLRASNAMMNNNDLVRKRRRLADLTGPIIPKCRSGV
Research Area	others
Source	E.coli
Target Names	TMOD1
Protein Names	Recommended name: Tropomodulin-1 Alternative name(s): Erythrocyte tropomodulin Short name= E-Tmod
Expression Region	1-359aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	67.6kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.