




# Recombinant Human Guanine nucleotide-binding protein subunit alpha-15 (GNA15)

<b>Product Code</b>	CSB-EP009587HU
<b>Relevance</b>	Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems.
<b>Abbreviation</b>	Recombinant Human GNA15 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P30679
<b>Alias</b>	Epididymis tissue protein Li 17E Guanine nucleotide-binding protein subunit alpha-16
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MARSLTWRCCPWCLTEDEKAAARVDQEINRILLEQKKQDRGELKLLLLGPGES GKSTFIKQMRIIHGAGYSEEERKGFRLVYQNIQVSMRAMIEAMERLQIPFSRP ESKHASLVMSQDPYKVTTFEKRYAAAMQWLWRDAGIRAYYERRREFHLLDS AVYYLSHLERITEEGYVPTAQDVLRSRMPPTTGINEYCFVQKTNLRIVDVGGQK SERKKWIHCFENVIALIYLASLSEYDQCLEENNQENRMKESLALFGTILELPWF KSTSVILFLNKTDILEEKIPTSHLATYFSPFQGPQDAEAAKRFILDMYTRMYTG CVDGPEGSKKGGARSRLFSHYTCATDTQNIRKVFKDVRDSVLARYLDEINLL
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	GNA15
<b>Protein Names</b>	Recommended name: Guanine nucleotide-binding protein subunit alpha-15 Short name= G alpha-15 Short name= G-protein subunit alpha-15 Alternative name(s): Guanine nucleotide-binding protein subunit alpha-16 Short name= G alpha-1
<b>Expression Region</b>	1-374aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-Trx-tagged
<b>Mol. Weight</b>	60.6kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.