




Recombinant Human F-actin-capping protein subunit beta (CAPZB)

Product Code	CSB-EP004514HU
Relevance	F-actin-capping proteins bind in a Ca ²⁺ -independent manner to the fast growing ends of actin filaments (barbed end) thereby blocking the exchange of subunits at these ends. Unlike other capping proteins (such as gelsolin and severin), these proteins do not sever actin filaments. Plays a role in the regulation of cell morphology and cytoskeletal organization.
Abbreviation	Recombinant Human CAPZB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47756
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSDQQLDCALDLMRRLPPQQIEKNLSLIDLVPCLCEDLLSSVDQPLKIARDKV VGKDYLLCDYNRDGDSYRSPWSNKYDPPLEDGAMP SARLRKLEVEANNAFD QYRDLYFEGGVSSVYLWDLDHGFAGVILIKKAGDGSKKIKGCWDSIHVVEVQE KSSGRTAHYKLTSTVMLWLQTNKSGSGTMNLGGS LTRQMEKDETVSDCSPHI ANIGRLVEDMENKIRSTLNEIYFGKTKDIVNGLRSVQTFADKSKQEALKNLVE ALKRKQQC
Research Area	Signal Transduction
Source	E.coli
Target Names	CAPZB
Protein Names	Recommended name: F-actin-capping protein subunit beta Alternative name(s): CapZ beta
Expression Region	1-272aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	57.6kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.