




Recombinant Mouse Tumor necrosis factor receptor superfamily member 4 (Tnfrsf4), partial

Product Code	CSB-EP023981MO
Relevance	Receptor for TNFSF4/OX40L/GP34. Is a costimulatory molecule implicated in long-term T-cell immunity (By similarity).
Abbreviation	Recombinant Mouse Tnfrsf4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47741
Alias	OX40 antigen OX40L receptor CD_antigen: CD134
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VTARRLNCVKHTYPSGHKCCRECQPGHGMVSRCDHTRDTLCHPCETGFYNE AVNYDTCKQCTQCNHRSGSELKQNCTPTQDTVCRCPGTQPRQDSGYKLG DCVPCPPGHFSPGNNQACKPWTNCTLSGKQTRHPASDSLDAVCEDRSLLATL LWETQRPTFRPTTVQSTTVWPRTSELPSPPTLVTPGEP
Research Area	Others
Source	E.coli
Target Names	Tnfrsf4
Expression Region	20-211aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	37.3kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.