



Recombinant Human Retinaldehyde dehydrogenase 3 (ALDH1A3)

| | |
|--------------------------|---|
| Product Code | CSB-EP001567HU-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P47895 |
| Storage Buffer | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Homo sapiens (Human) |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Sequence | ATANGAVEN GQPDRKPPAL PRPIRNLEVK FTKIFINNEW HESKSGKKFA TCNPSTREQI CEVEEGDKPD VDKAVEAAQV AFQRGSPWRR LDALSRGRLL HQLADLVERD RATLAALETM DTGKPFLHAF FIDLEGCIRT LRYFAGWADK IQGKTIPTDD NVVCFTRHEP IGVCGAITPW NFPLLMLVWK LAPALCCGNT MVLKPAEQTP LTALYLGSLI KEAGFPPGVV NIVPGFGPTV GAAISSHPQI NKIAFTGSTE VGKLVKEAAS RSNLKRVTLE LGGKNPCIVC ADADLDLAVE CAHQGVFFNQ GQCCTAASRV FVEEQVYSEF VRRSVEYAKK RPVGDPFDVK TEQGPQIDQK QFDKILELIE SGKKEGAKLE CGGSAMEDKG LFIKPTVFSE VTDNMRIAKE EIFGPVQPIL KFKSIEEVK RANSTDYGLT AAVFTKNLDK ALKLASALES GTVWINCYNA LYAQAPFGGF KMSGNGRELG EYALAEYTEV KTVTIKLGDK NP |
| Source | E.coli |
| Target Names | Name:ALDH1A3Synonyms:ALDH6 |
| Protein Names | Recommended name: Aldehyde dehydrogenase family 1 member A3 EC=1.2.1.5Alternative name(s): Aldehyde dehydrogenase 6 Retinaldehyde dehydrogenase 3 Short name= RALDH-3 Short name= RaLDH3 |
| Expression Region | 2-512 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.