




Recombinant Human Protein Wnt-7b (Wnt7B)

Product Code	CSB-EP026142HU
Relevance	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters
Abbreviation	Recombinant Human WNT7B protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56706
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ALSSVVALGANIICNKIPGLAPRQRAICQSRPDAIIVIGEGAQMGINECQYQFRF GRWNC SALGEKTVFGQELRVGSREAAFTYAITAAGVAHAVTAACSQGNLSNC GCDREKQGYYNQAEGWKWGGCSADVRYGIDFSRRFVDAREIKKNARRLMNL HNNEAGRKVLDRMQLECKCHGVSGSCTTKTCWTTLPKFREVGHLLKEKYNA AVQVEVVRASRLRQPTFLRIKQLRSYQKPMETDLVYIEKSPNYCEEDAATGSV GTQGRLCNRTSPGADGCDTMCCGRGYNTHQYTKVWQCNCCKFWCCFVKC NTCSERTEVFTCK
Research Area	Stem Cells
Source	E.coli
Target Names	WNT7B
Protein Names	Recommended name: Protein Wnt-7b
Expression Region	25-349aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Mol. Weight	56.3kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.