




Recombinant Human LIM domain transcription factor LMO4 (LMO4)

Product Code	CSB-EP013011HU
Relevance	Probable transcriptional factor.
Abbreviation	Recombinant Human LMO4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61968
Alias	Breast tumor autoantigen LIM domain only protein 4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MVNP GSSSQPPPVTAGSLSWKRCAGCGGKIADRFLLYAMDSYWHSRCLKCS CCQAQLGDIGTSCYTKSGMILCRNDYIRLFGNSGACSGQSIPASELVMRAQ GNVYHLKCFTCTCRNRLVPGDRFHYINGSLFCEHDRPTALINGHLNSLQSNP LLPDQKVC
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	LMO4
Protein Names	Recommended name: LIM domain transcription factor LMO4 Alternative name(s): Breast tumor autoantigen LIM domain only protein 4 Short name= LMO-4
Expression Region	1-165aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	45.0kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.