




Recombinant Human ADP-ribosylation factor 3 (ARF3)

Product Code	CSB-EP001990HU
Relevance	GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.
Abbreviation	Recombinant Human ARF3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61204
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GNIFGNLLKSLIGKKEMRILMVGLDAAGKTTILYKCLKLGEIVTTIPTIGFNVETVEY KNISFTVWDVGGQDKIRPLWRHYFQNTQGLIFVVDSDNRERVNEAREELMRM LAEDELRDVLLVFANKQDLPNAMNAAEITDKLGLHSLRHRNWYIQATCATSG DGLYEGLDWLANQLKNKK
Research Area	Signal Transduction
Source	E.coli
Target Names	ARF3
Protein Names	Recommended name: ADP-ribosylation factor 3
Expression Region	1-181aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	47.5kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.