




Recombinant Human Macrophage colony-stimulating factor 1 receptor (CSF1R), partial

Product Code	CSB-MP006044HU
Abbreviation	Recombinant Human CSF1R protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07333
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	IPVIEPSVPELVVKPGATVTLRCVGNNGSVEWDGPPSPHWTLYSDGSSSILSTN NATFQNTGTYRCTEPGDPLGGSAAIHLYVKDPAWPWNVLAQEVVVFEDQDAL LPCLLTDPVLEAGVSLVRVRGRPLMRHTNYSFSPWHGFTIHRAKFIQSQDYQC SALMGGRKVMSSIRLKVQKVIPGPPALTLVPAELVRIRGEAAQIVCSASSVDVN FDVFLQHNNTKLAIPQQSDFHNNRYQKVLTLNLDQVDFQHAGNYSVASNVQ GKHSTSMFFRVVESAYLNLSSQNLIQEVTVGEGNLKVMVEAYPGLQGFW TYLGPFSDHQPEPKLANATTKDTYRHTFTLSLPLKLPSEAGRYSFLARNPGGW RALTFELTLRYPPEVSVIWFINGSGTLLCAASGYQPQNVTLQCSGHTDRCD EAQVLQVWDDPYPEVLSQEPFHKVTVQSLTETLEHNQTYECRAHNSVGSG SWAFIPISAGAHTHPDE
Research Area	Cardiovascular
Source	Mammalian cell
Target Names	CSF1R
Protein Names	Recommended name: Macrophage colony-stimulating factor 1 receptor Alternative name(s): CSF-1 receptor Short name= CSF-1-R Short name= CSF-1R Short name= M-CSF-R EC= 2.7.10.1 Proto-oncogene c-Fms CD_antigen= C
Expression Region	20-512aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	82.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.