




# Recombinant Human Myosin light chain 3 (MYL3)

<b>Product Code</b>	CSB-EP015310HU
<b>Relevance</b>	Regulatory light chain of myosin. Does not bind calcium.
<b>Abbreviation</b>	Recombinant Human MYL3 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P08590
<b>Alias</b>	Cardiac myosin light chain 1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MAPKKPEPKKDDAKAAPKAAPAPAPPPEPERPKEVEFDASKIKIEFTPEQIEEF KEAFMLFDRTPKCEMKITYGQCGDVLRALGQNPTQAEVLRVLGKPRQEELNT KMMDFETFLPMLQHISKNKDTGTIEDFVEGLRVFDKEGNGTVMGAELRHVLA TLGERLTEDEVEKLMAGQEDSNGCINYEAFVKHIMSS
<b>Source</b>	E.coli
<b>Target Names</b>	MYL3
<b>Protein Names</b>	Recommended name: Myosin light chain 3 Alternative name(s): Cardiac myosin light chain 1 Short name= CMLC1 Myosin light chain 1, slow-twitch muscle B/ventricular isoform Short name= MLC1SB Ventricular/slow twitch myosin alkali light
<b>Expression Region</b>	1-195aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	48.8kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.