




Recombinant Human Defensin beta 4A (DEFB4A)

Product Code	CSB-EP006705HU(A4)
Relevance	Exhibits antimicrobial activity against Gram-negative bacteria and Gram-positive bacteria. May act as a ligand for C-C chemokine receptor CCR6. Can bind to both human and mouse CCR6 and induce chemotactic activity of CCR6-expressing cells
Abbreviation	Recombinant Human DEFB4A protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O15263
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MRVLYLLFSFLFIFLMPLPGVFVGGIGDPVTCLKSGAICHPVFCPRRYKQIGTCGLPGTKC CKKP
Source	E.coli
Target Names	DEFB4A
Protein Names	Recommended name: Beta-defensin 4AAlternative name(s): Beta-defensin 2 Short name= BD-2 Short name= hBD-2 Defensin, beta 2 Skin-antimicrobial peptide 1 Short name= SAP1
Expression Region	1-64aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	31.3kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.