



Recombinant Mouse Ubiquitin-protein ligase E3A (Ube3a), partial

Product Code	CSB-BP025488MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08759
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	NPADLKKQLYVEFEGEQGVDEGGVSKEFFQLVVEEIFNPDIGMFTYDEATKLF WFNPSSFETEGQFTLIGIVLGLAIYNNCILDVHFPMVVYRKLKMGKKGTFRDLDG SHPVLYQSLKDLLEYEGSVEDDMMITFQISQTDLFGNPMMYDLKENGDKIPITN ENRKEFVNLYSDYILNKSVEKQFKAFRRGFHMVTNESPLKYLFRPEEIELLIG SRNLDFQALEETTEYDGGYTRESVVIREFWEIFVHSFTDEQKRLFLQFTTGTDR APVGGGLGKLKMIIAKNGPDTERLPTSHTCFNVLLLPEYSSKEKLLKERLLKAITYA KGFGL
Source	Baculovirus
Target Names	Ube3a
Protein Names	Recommended name: Ubiquitin-protein ligase E3A EC= 6.3.2.-Alternative name(s): Oncogenic protein-associated protein E6-AP
Expression Region	542-870aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.