




# Recombinant Human Glutamyl-tRNA (Gln) amidotransferase subunit C, mitochondrial (GATC)

<b>Product Code</b>	CSB-EP009283HU
<b>Relevance</b>	Allows the formation of correctly charged Gln-tRNA(Gln) through the transamidation of misacylated Glu-tRNA(Gln) in the mitochondria. The reaction takes place in the presence of glutamine and ATP through an activated gamma-phospho-Glu-tRNA(Gln).
<b>Abbreviation</b>	Recombinant Human GATC protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O43716
<b>Alias</b>	Protein 15E1.2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MWSRLVWLGLRAPLGGRRQGFTSKADPQGSGRITAAVIEHLERLALVDFGSRE AVARLEKAIADFADRLRAVDTDGVEPMESVLEDRCLYLRSNDNVEGNCADPELLQ NSHRVVEEYFVAPPGNISLPLKDEQEPFPHS
<b>Research Area</b>	others
<b>Source</b>	E.coli
<b>Target Names</b>	GATC
<b>Protein Names</b>	Recommended name: Glutamyl-tRNA(Gln) amidotransferase subunit C, mitochondrial Short name= Glu-AdT subunit C EC= 6.3.5.-Alternative name(s): Protein 15E1.2
<b>Expression Region</b>	1-136aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	42.1kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.