



# Recombinant Human BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like (BNIP3L)

<b>Product Code</b>	CSB-BP002767HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O60238
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MSSHLVEPPPPLHNNNNNCEENEQSLPPPAGLNSSWVELPMNSSNGNDNGN GKNGGLEHVPSSSSSIHNGDMEKILLDAQHESGQSSSRGSSSHCDSPSPQEDG QIMFDVEMHTSRDHSSQSEEEVVEGEKEVEALKKSADWVSDWSSRPENIPPK EFHFRHPKRSVLSMRKSGAMKKGIFSAEFLKVFIPSLFLSHVLALGLGIYIGK RLSTPSASTY
<b>Source</b>	Baculovirus
<b>Target Names</b>	BNIP3L
<b>Protein Names</b>	Recommended name: BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like Alternative name(s): Adenovirus E1B19K-binding protein B5 BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A NIP3-like protein X Short name= NIP3L
<b>Expression Region</b>	1-219
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.