



Recombinant Human BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like (BNIP3L)

| | |
|--------------------------|---|
| Product Code | CSB-EP002767HU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | O60238 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Sequence | MSSHLVEPPPPLHNNNNNCEENEQSLPPPAGLNSSWVELPMNSSNGNDNGN GKNGGLEHVPSSSSSIHNGDMEKILLDAQHESGQSSSRGSSSHCDSPSPQEDG QIMFDVEMHTSRDHSSQSEEEVVEGEKEVEALKKSADWVSDWSSRPENIPPK EFHFRHPKRSVLSMRKSGAMKKGIFSAEFLKVFIPSLFLSHVLALGLGIYIGK RLSTPSASTY |
| Source | E.coli |
| Target Names | BNIP3L |
| Protein Names | Recommended name: BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like Alternative name(s): Adenovirus E1B19K-binding protein B5 BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A NIP3-like protein X Short name= NIP3L |
| Expression Region | 1-219 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |