




Recombinant Human Dickkopf-related protein 1 (DKK1), partial Escherichia Coli

Product Code	CSB-EP006920HU
Abbreviation	Recombinant Human DKK1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94907
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	LNSVLNSN AIKNLPPPLG GAAGHPGSAV SAAPGILYPG GNKYQTIDNY QPYPCAEED E CGTDEYCASP TRGGDAGVQI CLACRKRKR CMRHAMCCPG NYCKNGICVS SDQNHFRGEI EETITESFGN DHSTLDGYSR RTTLSSKMYH TKGQEGSVCL RSSDCASGLC CARHFWSKIC KPVLKEGQVC TKHRRKGSYG LEIFQRCYCG EGLSRIQKD HHQASNSSRL HTCQRH
Source	E.coli
Target Names	DKK1
Protein Names	Recommended name: Dickkopf-related protein 1 Short name= Dickkopf-1 Short name= Dkk-1 Short name= hDkk-1Alternative name(s): SK
Expression Region	33-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Protein Length	partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.