



# Recombinant Rat Vitamin K-dependent protein C (Proc)

<b>Product Code</b>	CSB-BP018742RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P31394
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ANSFLEEVR AGSLERECME EICDFEEAQE IFQNVEDTLA FWIKYFDGDQ CSTPPLDHQC DSPCCGHGTC IDGLGGFSCS CDKGWEGRFC QQEMGFQDCR VKNGGCYHYC LEETRGRRCR CAPGYELADD HMHCRPTVNF PCGKLWKRTD KKRKNFKRDI DPEDEELELG PRIVNGTLTK QGDSPWQAIL LDSKKKLACG GVLHTSWVL TAAHCLESSK KLTVRLGEYD LRRRDPWELD LDIKEVLVHP NYTRSNSDND IALLRLSQPA TLSKTIVPIC LPNSGLAQEL SQAGQETVVT GWGYQSDKVK DGRRNRTFIL TFIRIPLAAR NDCMQVMNNV VSENMLCAGI IGDTRDACDG DSGGPMVFFF RGTWFLVGLV SWGEGCGHLN NYGVYTKVGS YLKWIIHSYIG ERDVSLKSPK L
<b>Source</b>	Baculovirus
<b>Target Names</b>	Proc
<b>Protein Names</b>	Recommended name: Vitamin K-dependent protein C EC= 3.4.21.69 Alternative name(s): Anticoagulant protein C Autoprothrombin IIA Blood coagulation factor XIV Cleaved into the following 3 chains: 1. Vitamin K-dependent protein C li
<b>Expression Region</b>	42-461
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.