




Recombinant Human Histone H1.0 (H1-0)

Product Code	CSB-EP010087HU
Relevance	Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The H1F0 histones are found in cells that are in terminal stages of differentiation or that have low rates of cell division.
Abbreviation	Recombinant Human H1-0 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07305
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MTENSTSAPAAKPKRAKASKKSTDHPKYSDMIVAAIQAEKNRAGSSRQSIQKYI KSHYKVGGENADSQIKLSIKRLVTTGVLKQTKGVBGASGSFRLAKSDEPKKSVAF KKTKEIKKVATPKKASKPKKAASKAPTCKPKATPVKKAKKKLAATPKKAKKPK TVKAKPVKASKPKKAKPKKAKSSAKRAGKKK
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	H1-0
Protein Names	Recommended name: Histone H1.0Alternative name(s): Histone H1' Histone H1(0)
Expression Region	1-194aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-GST-tagged
Mol. Weight	52.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.