




Recombinant Pan troglodytes Acid ceramidase (ASAH1)

Product Code	CSB-EP002169EQV
Relevance	Hydrolyzes the sphingolipid ceramide into sphingosine and free fatty acid.
Abbreviation	Recombinant Chimpanzee ASAH1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5A6P2
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QHAPPWTEDCRKSTYPPSGPTYRGPVWPWYTINLDLPPYKRWHELMMLDKAPML KVIVNSLKNMINTFVPSGKIVQVVDEKLPGLLGNFPGPFFEEEMKGIAAVTDIPLG EIIISFNIFYELFTICTSIVAEDKKGHLIHGRNMDFGVFLGWNINNDTWVITEQLKP LTVNLDFQRNNKTVFKASSFAGYVGMMLTGFKPGLFSLSLNERFSINGGYLGILE WILGKKDAMWIGFLTRTVLENSTSYEEAKNLLTKTKILAPAYFILGGNQSSEGC VITRDRKESLDVYELDAKQGRWYVQTNYDRWKHPFFLDDRRTPAKMCLNRT SQENISFETMYDVLSTKPVLNKLTVYTTLIDVTKGQFETYLRDCPDPCIGW
Research Area	Signal Transduction
Source	E.coli
Target Names	ASAH1
Protein Names	Acylsphingosine deacylase N-acylsphingosine amidohydrolase
Expression Region	22-395aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	47.7 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.