




Recombinant Mycoplasma pneumoniae ADP-ribosylating toxin CARDS (cards), partial

Product Code	CSB-EP300476MLW
Abbreviation	Recombinant Mycoplasma pneumoniae cards protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P75409
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MPNPVRFVYRVDLRSPEEIFEHGFSTLGDVRNFFEHILSTNFGRSYFISTSETP TAAIRFFGSWLREYVPEHPRRAYLYEIRADQHFYNARATGENLLDLMRQRQVV FDSGDREMAQM GIRALRTSFAYQREWFTDGPAAANVRS AWLVDAVPVEPGH AHPAGRVVETTRINEPEMHNP HYQELQTQANDQPWLPTPGIA
Research Area	Others
Source	E.coli
Target Names	cards
Protein Names	Recommended name: ADP-ribosylating toxin CARDS EC= 2.4.2.-Alternative name(s): ADP-ribosyltransferase CARDS CARDX TX
Expression Region	1-202aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	30.4 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.