



Recombinant Human ETS translocation variant 4 (ETV4)

Product Code	CSB-BP007858HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43268
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MERRMKAGYL DQQVPYTFSS KSPGNGLRE ALIGPLGKLM DPGSLPPLDS EDLFQDLSHF QETWLAEAQV PDSDEQFVPD FHSENLAHFS PTTRIKKEPQ SPRTDPALSC SRKPPLPYHH GEQCLYSSAY DPPRQIAIKS PAPGALGQSP LQPFPRAEQR NFLRSSGTSQ PHPGHGYLGE HSSVFQQPLD ICHSFTSQGG GREPLPAPYQ HQLSEPCPPY PQQSFKQEYH DPLYEQAGQP AVDQGGVNGH RYPGAGVVIK QEQTDFAYDS DVTGCASMYL HTEGFSGPSP GDGAMGYGYE KPLRPFDDV CVVPEKFEGD IKQEGVGAFR EGPPYQRRGA LQLWQFLVAL LDDPTNAHFI AWTGRGMEFK LIEPEEVARL WGIQKNRPAM NYDKLSRSLR YYYEKGIMQK VAGERYVYKF VCEPEALFSL AFPDNQRPAL KAEDRPVSE EDTVPLSHLD ESPAYLPELA GPAQPFPGPKG GYSY
Source	Baculovirus
Target Names	ETV4
Protein Names	Recommended name: ETS translocation variant 4 Alternative name(s): Adenovirus E1A enhancer-binding protein E1A-F Polyomavirus enhancer activator 3 homolog Short name= Protein PEA3
Expression Region	1-484
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.