



# Recombinant Human ETS translocation variant 4 (ETV4)

<b>Product Code</b>	CSB-YP007858HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P43268
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MERRMKAGYL DQQVPYTFSS KSPGNGSLRE ALIGPLGKLM DPGSLPPLDS EDLFQDLSHF QETWLAEAQV PDSDEQFVPD FHSENLAHFS PTTRIKKEPQ SPRTDPALSC SRKPPLPYHH GEQCLYSSAY DPPRQIAIKS PAPGALGQSP LQPFPRAEQR NFLRSSGTSQ PHPGHGYLGE HSSVFQQPLD ICHSFTSQGG GREPLPAPYQ HQLSEPCPPY PQQSFKQEYH DPLYEQAGQP AVDQGGVNGH RYPGAGVVIK QEQTDFAYDS DVTGCASMYL HTEGFSGPSP GDGAMGYGYE KPLRPFDDV CUVPEKFEGD IKQEGVGAFR EGPPYQRRGA LQLWQFLVAL LDDPTNAHFI AWTGRGMEFK LIEPEEVARL WGIQKNRPAM NYDKLSRSLR YYYEKGIMQK VAGERYVYKF VCEPEALFSL AFPDNQRPAL KAEDRPVSE EDTVPLSHLD ESPAYLPELA GPAQPFPGPKG GYSY
<b>Source</b>	Yeast
<b>Target Names</b>	ETV4
<b>Protein Names</b>	Recommended name: ETS translocation variant 4 Alternative name(s): Adenovirus E1A enhancer-binding protein E1A-F Polyomavirus enhancer activator 3 homolog Short name= Protein PEA3
<b>Expression Region</b>	1-484
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.