



Recombinant Mouse Histone H1.4 (H1-4)

Product Code	CSB-EP010380MO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43274
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SETAPAAPA APAPAEKTPV KKKARKAAGG AKRKTSGPPV SELITKAVAA SKERSGVSLA ALKKALAAAG YDVEKNNSRI KLGLKSLVSK GTLVQTKGTG ASGSFKLNKK AASGEAKPKA KRAGAACKAKK PAGAACKPKK AAGTATAKKS TKKTPKKAKK PAAAAGAKKA KSPKKAKATK AKKAPKSPAK AKTVKPKAAK PKTSKPKAAK PKKTAAKKK
Source	E.coli
Target Names	H1-4
Protein Names	Recommended name: Histone H1.4Alternative name(s): H1 VAR.2 H1e
Expression Region	2-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.