




Recombinant Avena sativa Endochitinase, partial

Product Code	CSB-BP308852DPO
Relevance	This protein functions as a defense against chitin-containing fungal pathogens.
Abbreviation	Recombinant Avena sativa Endochitinase protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P86181
Product Type	Recombinant Proteins
Immunogen Species	Avena sativa (Oat)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VSSVISSSLFEKMLLHRGFYTYDAFIAAAKSFPAPATTGSTDVRKREVA AFLAQ TSHETTGGWPTAPDGPYELGSTSDYFGRGPIQISYNNYGAAGKAIGVDLLRN PDLVTS DNTVEFKTALWFWMT PQSPKPSSHDVITGRWSPSSTDKAAGRVPGY GVL TNIIDGGVECGKGQESHVADRIGYYKDNLDCYNQKPFA
Research Area	others
Source	Baculovirus
Target Names	N/A
Protein Names	Recommended name: Endochitinase EC= 3.2.1.14
Expression Region	1-200aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	25.7kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.