




Recombinant Dactylis glomeRata Pollen allergen Dac g 3, partial

Product Code	CSB-EP309006DAC
Abbreviation	Recombinant Dactylis glomerata Pollen allergen Dac g 3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P93124
Product Type	Recombinant Proteins
Immunogen Species	Dactylis glomerata (Orchard grass) (Cock's-foot grass)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VKVTFKVEKGSDPKKLVLDIKYTRPGDTLAEVELRQHGSEEWELTKKGNLWE VKSSKPLTGPFNFRFMSKGGMRNVFDEVIPTAFKIGTTYTPEE
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Pollen allergen Dac g 3 Alternative name(s): Allergen Dac g III Allergen= Dac g 3
Expression Region	1-96aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	26.9kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.