



Recombinant Human Melanoma-associated antigen 11 (MAGEA11)

Product Code	CSB-MP013325HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43364
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	METQFRRGGL GCSPASIKRK KKREDSGDFG LQVSTMFSED DFQSTERAPY GPQLQWSQDL PRVQVFREQA NLEDSPRRT QRITGGEQVL WGPITQIFPT VRPADLTRVI MPLEQRSQHC KPEEGLQAQE EDLGLVGAQA LQAEEQEAAF FSSTLNVGTL EELPAAESPS PPQSPQEESF SPTAMDAIFG SLSDEGSGSQ EKEGPSTSPD LIDPESFSQD ILHDKIIDLV HLLLRKYRVK GLITKAEMLG SVIKNYEDYF PEIFREASVC MQLLFGIDVK EVDPTSHSYV LVTSLNLSYD GIQCNEQSMP KSGLLIIVLG VIFMEGNCIP EEVMWEVLSI MGVYAGREHF LFGPEKRLLT QNWVQEKEYLV YRQVPGTDP A CYEFLWGPRA HAETSKMKVL EYIANANGRD PTSYPSLYED ALREEGEGV
Source	Mammalian cell
Target Names	MAGEA11
Protein Names	Recommended name: Melanoma-associated antigen 11 Alternative name(s): Cancer/testis antigen 1.11 Short name= CT1.11 MAGE-11 antigen
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.