



Recombinant Abrus precatorius Abrin-a, partial

Product Code	CSB-BP319970AAC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11140
Product Type	Recombinant Protein
Immunogen Species	Abrus precatorius (Indian licorice) (Glycine abrus)
Purity	>85% (SDS-PAGE)
Sequence	QDRPIKFSTE GATSQSYKQF IEALRERLRG GLIHDIPVLP DPTTLQERNR YITVELSNSD TESIEVGIDV TNAYVVAYRA GTQSYFLRDA PSSASDYLFT GTDQHSLPFY GTYGDLERWA HQSRQQIPLG LQALTHGISF FRSGGNDNEE KARTLIVIIQ MVAEAARFRY ISNRVRVSIQ TGTAFAQPDAA MISLENNWDN LSRGVQESVQ DTFPNQVTLT NIRNEPVIVD SLSHPTVAVL ALMLFVCNPP N
Source	Baculovirus
Protein Names	Recommended name: Abrin-a Cleaved into the following 3 chains: 1. Abrin-a A chain EC= 2. 3.2.2.22 Alternative name(s): rRNA N-glycosidase Linker peptide Abrin-a B chain
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.