




# Recombinant Neisseria meningitidis serogroup B Major outer membrane protein P.IA (porA)

<b>Product Code</b>	CSB-EP320687NGG
<b>Relevance</b>	Serves as a slightly cation selective porin. Major antigen on the gonococcal cell surface and it may have pathogenic properties in addition to its porin activity.
<b>Abbreviation</b>	Recombinant Neisseria meningitidis serogroup B porA protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0DH58
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Neisseria meningitidis serogroup B (strain MC58)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	DVSLYGEIKAGVEGRNYQLQLTEAQAANGGASGQVKVTKVTKAKSRIRTKISD FGSFIGFKGSEDLGDGLKAVWQLEQDVSVAGGGATQWGNRESFIGLAGEFGT LRAGR VANQFDDASQAIDPWDSNNDVASQLGIFKRHDDMPVSVRYDSPEFSG FSGSVQFVPIQNSKSAYTPAYYTKNTNNLTLVPAVVGKPGSDVYYAGLNYKN GGFAGNYAFKYARHANVGRNAFELFLIGSGSDQAKGTDPLKNHQVHRLTGGY EEGGLNLALAAQLDLENGDKTKNSTTEIAATASYRFGNAVPRISYAHGFDFIE RGKKGENTS YDQIIAGVDYDFSKRTSAIVSGAWLKRNTGIGNYQTQINAASVGLR HKF
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	porA
<b>Protein Names</b>	Class 1 protein
<b>Expression Region</b>	20-392aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	44.1 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.