



Recombinant Sheep Dipeptidase 1 (DPEP1)

Product Code	CSB-BP007123SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43477
Product Type	Recombinant Protein
Immunogen Species	Ovis aries (Sheep)
Purity	≥85% (SDS-PAGE)
Sequence	DQFR DNAVRLMQST PVIDGHNDLP WALLKKFNNQ LQDPRANLTS LNSTHTNIPK LKAGFVGAQF WSAYTPCDTQ NKDSVKRTVE QIDVIQRMCC LYPETFLCVT DSAGIQQAFQ EGKVASLVGV EGGHSIDSSL GVLRALYHLG MRYLTLTHSC NTPWADNWL V DTGEDKAQSQ GLSSFGQSVV KEMNRLGIII DLAHVSVATM EAALQLSKAP VIFSHSSAYS LCHHRRNVPD HVLQLVKQTG SLVMVNFYND YVSCAEANL SQVADHLDYI KKVAGAGAVG FGGDYDGVSR LPSGLEDVSK YPDLVAELLR RQWTEEEVRG ALAENLLRVF KAVEQASDHK QAPGEEPIPL GQLEASCRTK YGYS
Source	Baculovirus
Target Names	DPEP1
Protein Names	Recommended name: Dipeptidase 1 EC= 3.4.13.19 Alternative name(s): Microsomal dipeptidase
Expression Region	17-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.