






Recombinant Human HLA class I histocompatibility antigen, A alpha chain (HLA-A), partial

Product Code	CSB-EP328989HU
Relevance	Involved in the presentation of foreign antigens to the immune system.
Abbreviation	Recombinant Human HLA-A protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SUL5
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	GSHSMRYFFTSVSRPGRGEPFRFIAVG YVDDTQFVRFDSDAASQKMEPRAPWI EQEGPEYWDQETRNMKAHSQTD RANLGLRGYYNQSE DGSHTIQIMYGCDV GPDGRFLRGYRQDAYDGKD YIALNEDLRSWTAADMAAQITKRKWEAVHAAE QRRVYLEGRCVDGLRRYLENGKETLQRTDPPKTHMTHHPISDHEATLRCWAL GFYPAEITLTWQRDGEDQTQDTELVETRPAGDGT FQKWAAVVVP SGEEQRY TCHVQHEGLPKPLTLRWELSSQPTIPI
Research Area	Others
Source	E.coli
Target Names	HLA-A
Protein Names	MHC class I antigen A*1HLAA
Expression Region	25-308aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	36.7 kDa
Protein Length	Extracellular Domain
Image	



The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.