




Recombinant Naja mossambica Cytotoxin 5

Product Code	CSB-EP329557NAH
Relevance	Shows cytolytic activity on many different cells by forming pore in lipid membranes. In vivo, increases heart rate or kills the animal by cardiac arrest. In addition, it binds to heparin with high affinity, interacts with Kv channel-interacting protein 1 (KCNIP1) in a calcium-independent manner, and binds to integrin alpha-V/beta-3 (ITGAV/ITGB3) with moderate affinity.
Abbreviation	Recombinant Naja mossambica Cytotoxin 5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25517
Product Type	Recombinant Protein
Immunogen Species	Naja mossambica (Mozambique spitting cobra)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	LKCKKKLIPLFSKTCPEGKNLCYKMTMRLAPKVPVKRGCIDVCPKSSFLVKYECCDTDRCN
Research Area	others
Source	E.coli
Protein Names	CTX M51 CTX V
Expression Region	1-60aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	10.8 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.