




# Recombinant Crotalus adamanteus Snake venom metalloproteinase adamalysin-2

<b>Product Code</b>	CSB-EP330325DYB
<b>Relevance</b>	Has no significant hemorrhagic activity, but inactivates serpins by limited proteolysis of their reactive-site loops.
<b>Abbreviation</b>	Recombinant Crotalus adamanteus Snake venom metalloproteinase adamalysin-2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34179
<b>Alias</b>	Adamalysin II Proteinase II
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Crotalus adamanteus (Eastern diamondback rattlesnake)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	QQNLPQRYIELVVADRRVFMKYNSDLNIIRTRVHEIVNIINGFYRSLNIDVSLVN LEIWSGQDPLTIQSSSNTLNSEGLWREKVLLNKKKKDNAQLLTAIEFKCETLG KAYLNSMCMNRSSVGIVKDHSPINLLVAVTMAHELGHNLGMEHDGKDCLRGA SLCIMRPGLTPGRSYEFSDDSMGYYQKFLNQYKPQCILNKP
<b>Source</b>	E.coli
<b>Protein Names</b>	Recommended name: Zinc metalloproteinase adamalysin-2 EC=3.4.24.46 Alternative name(s): Adamalysin II Proteinase II
<b>Expression Region</b>	1-203aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	27.1kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.