



Recombinant Clostridium botulinum Hemagglutinin component of the neurotoxin complex (ha17)

Product Code	CSB-MP332826CLQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5HZZ5
Product Type	Recombinant Protein
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SSERTFLPN GNYKIKSLFS DSLYLTYSYG ALSFSNTSSL DNQKWKLEYI SSSNGFRFSN VAEPNKYLAY NDYGFYLLSS SSNNSLWNPI KIAINSYIIC TLSIVNVTDY AWTIYDNNNN ITDQPILNLP NFDINNSNQI LKLEKL
Source	Mammalian cell
Protein Names	Recommended name: Hemagglutinin component HA-17 Alternative name(s): HA 17 kDa subunit
Expression Region	2-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.