



Recombinant Vespa vulgaris Phospholipase A1

Product Code	CSB-BP343003VET
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49369
Product Type	Recombinant Protein
Immunogen Species	Vespa vulgaris (Yellow jacket) (Wasp)
Purity	>85% (SDS-PAGE)
Sequence	GPKC PFNSD TVSII IETREN RNRD LYTLQ TLQNH PEFK KKTITR PVVFITHGFT SSASE TNFIN LAKALVDKDN YMVISIDWQT AACTNEAAGL KYL YPTAAR NTRLV GQYIA TITQKLVKHY KISMANIRLI GHSLGAHASG FAGKKVQELK LGKYSEIIGL DPARPSFDSN HCSERLCETD AEYVQIIHTS NYLGTEKTLG TVDFYM NNGK NQPGCGRFFS EVCSH SRAVI YMAECIKHEC CLIGIPKSKS SQPISSCTKQ ECVCVGLNAK KYPSRGSFYV PVESTAPFCN NKGKII
Source	Baculovirus
Protein Names	Recommended name: Phospholipase A1 EC= 3.1.1.32 EC= 3.1.1.4 Alternative name(s): Allergen Ves v I Allergen= Ves v 1
Expression Region	37-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.