



Recombinant Centruroides noxius Beta-mammal toxin Cn2

Product Code	CSB-MP355662DRB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01495
Product Type	Recombinant Protein
Immunogen Species	Centruroides noxius (Mexican scorpion)
Purity	>85% (SDS-PAGE)
Sequence	KEGY LVDKNTGCKY ECLKLGDNDY CLRECKQQYG KGAGGYCYAF ACWCTHLYEQ AIVWPLPNKR CS
Source	Mammalian cell
Protein Names	Recommended name: Beta-mammal toxin Cn2 Short name= Toxin 2 Alternative name(s): Toxin II.9.2.2
Expression Region	17-82
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.