




Recombinant Triticum aestivum Glutenin, high molecular weight subunit PC237, partial

Product Code	CSB-EP355887TQN
Relevance	Glutenins are high-molecular weight seed storage proteins of wheat endosperm. Thought to be responsible for the visco-elastic property of wheat dough.
Abbreviation	Recombinant Triticum aestivum Glutenin, high molecular weight subunit PC237 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P02862
Product Type	Recombinant Proteins
Immunogen Species	Triticum aestivum (Wheat)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	LVSVEHQAARLKVAKAQQLAAQLPAMCRLEGGDALSASQ
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Glutenin, high molecular weight subunit PC237
Expression Region	1-39aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	31.3kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.