




Recombinant Epstein-Barr virus Envelope glycoprotein L (gL)

Product Code	CSB-EP355982EFA
Relevance	The heterodimer glycoprotein H-glycoprotein L is required for the fusion of viral and plasma membranes leading to virus entry into the host cell. Membrane fusion is mediated by the fusion machinery composed at least of gB and the heterodimer gH/gL. Fusion of EBV with B-lymphocytes requires the additional receptor-binding protein gp42, which forms a complex with gH/gL. May also be required for virus attachment to epithelial cells
Abbreviation	Recombinant Epstein-Barr virus gL protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03212
Alias	Glycoprotein 25 Short name: gp25
Product Type	Recombinant Protein
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	NWAYPCCHVTQLRAQHLLALENISDIYLVSNQTCDFSLASLNSPKNGSNQLV SRCANGLNVVVSFFISILKRSSSALTGHLRELLTTLETLYGFSFVEDLFGANLNRY AWHRGG
Research Area	Others
Source	E.coli
Target Names	gL
Expression Region	23-137aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	28.7kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.