




Recombinant Mouse Low affinity immunoglobulin gamma Fc region receptor III (Fcgr3), partial

Product Code	CSB-EP357412MO
Relevance	Receptor for the Fc region of complexed immunoglobulins gamma. Low affinity receptor.
Abbreviation	Recombinant Mouse Fcgr3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08508
Alias	Fc-gamma RIII Short name: FcRIII CD_antigen: CD16
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ALPKAVVKLDPPWVQLKEDMTLMCEGTHNPGNSSTQWFHNGRSIRSQVQA SYTFKATVNDSGEYRCQMEQTRLSDPVDLGVISDWLLLQTPQRFLEGETITL RCHSWRNKLLNRISFFHNEKSVRYHHYKSNFSIPKANHSHSGDYCKGSLGS TQHQSKPVTITVQDPATTSSISLVWYHT
Research Area	Others
Source	E.coli
Target Names	Fcgr3
Expression Region	31-215aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	37.2kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.