



# Recombinant Escherichia coli O6:H1 Exopolyphosphatase (ppx)

<b>Product Code</b>	CSB-EP360285EGX-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0AFL7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PIHDKSPRP QEFAAVDLGS NSFHMVIARV VDGAMQIIGR LKQRVHLADG LGPDNMLSEE AMTRGLNCLS LFAERLQGFS PASVCIVGTH TLRQALNATD FLKRAEKVIP YPIEISGNE EARLIFMGVE HTQPEKGRKL VIDIGGGSTE LVIGENFEPI LVESRRMGCV SFAQLYFPGG VINKENFQRA RMAAAQKLET LTWQFRIQGW NVAMGASGTI KAAHEVLMEM GEKDGIITPE RLEKLVKEVL RHRNFASLSL PGLSEERKTV FVPGLAILCG VFDALAIREL RLSGDALREG VLYEMEGFRFR HQDVRSRTAS SLANQYHIDS EQARRVLDTT MQMYEQWREQ QPKLAHPQLE ALLRWAAMLH EVGLNINHSG LHRHSAYILQ NSDLPGFNQE QQLMMATLVR YHRKAIKLLD LPRFTLFKKK QFLPLIQLLR LGVLLNNQRQ ATTTPTTLTL ITDDSHWTLR FPHDWFSQNA LVLLDLEKEQ EYWEGVAGWR LKIEEESTPE IAA
<b>Source</b>	E.coli
<b>Target Names</b>	ppx
<b>Protein Names</b>	Recommended name: Exopolyphosphatase Short name= ExopolyPase EC= 3.6.1.11Alternative name(s): Metaphosphatase
<b>Expression Region</b>	2-513
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.