




# Recombinant Enterobacter aerogenes Ribitol 2-dehydrogenase (rbtD)

<b>Product Code</b>	CSB-EP360400EKN
<b>Abbreviation</b>	Recombinant Enterobacter aerogenes rbtD protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P00335
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Klebsiella aerogenes (Enterobacter aerogenes)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MKHSVSSMNTSLSGKVAAITGAASGIGLECARTLLGAGAKVVLIDREGEKLNKL VAELGENAFALQVDLMQADQVDNLLQGILQLTGRLDIFHANAGAYIGGPVAEG DPDVWDRVLHLNINAAFRCVRSVLPHLIAQKSGDIIFTAVIAGVVPVIWEPVYTA SKFAVQAFVHTTRRQVAQYGVRVGAVLPGPVVTALLDDWPKAKMDEALANG SLMQPIEVAESVLFMVTRSKNVTVRDIVILPNSVDL
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	rbtD
<b>Protein Names</b>	Recommended name: Ribitol 2-dehydrogenase Short name= RDH EC= 1.1.1.56
<b>Expression Region</b>	1-249aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	30.5 kDa
<b>Protein Length</b>	Full Length

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.