




Recombinant Escherichia coli Tetracycline repressor protein class B from transposon Tn10 (tetR)

Product Code	CSB-EP366163ENL
Relevance	TetR is the repressor of the tetracycline resistance element; its N-terminal region forms a helix-turn-helix structure and binds DNA. Binding of tetracycline to TetR reduces the repressor affinity for the tetracycline resistance gene (tetA) promoter operator sites.
Abbreviation	Recombinant E.coli tetR protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04483
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSRLDKSKVINSALELLNEVGIEGLTTRKLAQKLGVEQPTLYWHVKNKRALLDA LAIEMLDLRHHTHFCPLEGESWQDFLRNNAKSFRCALLSHRDGAKVHLGTRPT EKQYETLENQLAFLCQQGFSLLENALYALSAVGHFTLGCVLEDQEHQVAKEER ETPTTDSMPPLLRQAIELFDHQGAEPFLFGLELIICGLEKQLKCESGS
Research Area	Others
Source	E.coli
Target Names	tetR
Protein Names	Recommended name: Tetracycline repressor protein class B from transposon Tn10
Expression Region	1-207aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	39.4kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.