



Recombinant Mouse Catalase (Cat)

Product Code	CSB-MP004563MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24270
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p>SDSRDPASD QMKQWKEQRA SQRPDVLTGG GGNPIGDKLN IMTAGSRGPL LVQDVVFTDE MAHFDREIRIP ERVVHAKGAG AFGYFEVTHD ITRYSKAKVF EHIGKRTPIA VRFSTVTGES GSADTVRDPG GFAVKFYTED GNWDLVGNNT PIFFIRDAIL FPSFIHSQKR NPQTHLKDPD MVWDFWVSLRP ESLHQVSFLF SDRGIPDGHR HMNGYGSHTF KLVNADGEAV YCKFHYKTDQ GIKNLPVGEA GRLAQEDPDY GLRDLFNAIA NGNYPSWTFY IQVMTFKEAE TFPFNPFDLT KVWPHKDYPL IPVGKLVLNK NPVNYFAEVE QMAFDPSNMP PGIEPSPDKM LQGRLFAYPD THRHRLGPNY LQIPVNCPYR ARVANYQRDG PMCMHDNQQG APNYYPNSFS APEQQRSALE HSVQCAVDVK RFNSANEDNV TQVRTFYTKV LNEEERKRLC ENIAGHLKDA QLFIQKKAVK NFTDVHPDYG ARIQALLDKY NAEKPKNAIH TYTQAGSHMA AKGKANL</p>
Source	Mammalian cell
Target Names	Cat
Protein Names	Recommended name: Catalase EC= 1.11.1.6
Expression Region	2-527
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	<p>This gene encodes catalase, a key antioxidant enzyme in the bodies defense against oxidative stress. Catalase is a heme enzyme that is present in the peroxisome of nearly all aerobic cells. Catalase converts the reactive oxygen species hydrogen peroxide to water and oxygen and thereby mitigates the toxic effects of hydrogen peroxide. Oxidative stress is hypothesized to play a role in the development of many chronic or late-onset diseases such as diabetes, asthma, Alzheimer s disease, systemic lupus erythematosus, rheumatoid arthritis, and cancers. Polymorphisms in this gene have been associated with decreases in catalase activity but, to date, acatalasemia is the only disease known to be caused by this gene.</p>
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final</p>



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.