



Recombinant Influenza A virus Matrix protein 1 (M)

Product Code	CSB-EP389592ILS
Relevance	Plays critical roles in virus replication, from virus entry and uncoating to assembly and budding of the virus particle. M1 binding to ribonucleocapsids (RNPs) in nucleus seems to inhibit viral transcription. Interaction of viral NEP with M1-RNP is thought to promote nuclear export of the complex, which is targeted to the virion assembly site at the apical plasma membrane in polarized epithelial cells. Interactions with NA and HA may bring M1, a non-raft-associated protein, into lipid rafts. Forms a continuous shell on the inner side of the lipid bilayer in virion, where it binds the RNP. During virus entry into cell, the M2 ion channel acidifies the internal virion core, inducing M1 dissociation from the RNP. M1-free RNPs are transported to the nucleus, where viral transcription and replication can take place Determines the virion's shape: spherical or filamentous. Clinical isolates of influenza are characterized by the presence of significant proportion of filamentous virions, whereas after multiple passage on eggs or cell culture, virions have only spherical morphology. Filamentous virions are thought to be important to infect neighboring cells, and spherical virions more suited to spread through aerosol between hosts organisms
Abbreviation	Recombinant Influenza A virus Matrix protein 1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A4GCLO
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/USA:Iowa/1943 H1N1)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSLLTEVETYVLSIVPSGPLKAEIAQRLEDVFAGKNTDLEALMEWLKTRPILSPL TKGILGFVFTLTVPSEGLQRRRFVQNALNGNGDPNNMDRAVKLYRKLKREIT FHGAKEIALSYSAGALASCMGLIYNRMGAVTTEVAFGLVCATCEQIADSQHRS HRQMVTNTNPLIRHENRMVLASTTAKAMEQMAGSSEQAAEAMEVASQARQM VQAMRAIGTHPSSSAGLKNLLENLQAYQKRMGVQMQRFK
Research Area	others
Source	E.coli
Target Names	M
Protein Names	Recommended name: Matrix protein 1 Short name= M1
Expression Region	1-252aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged

