



# Recombinant Mouse Radixin (Rdx)

<b>Product Code</b>	CSB-YP019532MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P26043
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>MPKPINVRVT TMDAELEFAI QPNTTGKQLF DQVVKTVGLR EVWFFGLQYV            DSKGYSTWLK LNKKVTQQDV KKENPLQFKF RAKFFPEDVS EELIQEITQR            LFFLQVKEAI LNDEIYCPPE TAVLLASYAV QAKYGDYNKE IHKPGYLAND            RLLPQRVLEQ HKLTKEQWEE RIQNWHEEHR GMLREDSMME YLKIAQDLEM            YGVNYFEIKN KKGTELWLGV DALGLNIYEH DDKLTPKIGF PWSEIRNISF            NDKKFKVIKPI DKKAPDFVfy APRLRINKRI LALCMGNHEL YMRRRKPDTI            EVQQMKAQAR EEKHQKQLER AQLENEKKKR EIAEKEKERI EREKEELMER            LRQIEEQTVK AQKELEEQTR KALELEQERQ RAKEEAERLD RERRAAEEAK            SAIKQAADQ MKNQEQLAAE LAEFTAKIAL LEEAKKKKKEE EATEWQHKAF            AAQEDLEKTK EELKTVMSAP PPPPPPPVIP PTENEHDEQD ENSAEASAE            SSEGVMNHRs EEERTVETQK NERVKKQLQA LSSELAQARD ETKKTQNDVL            HAENVKAGRD KYKTLRQIRQ GNTKQRIDEF EAM</p>
<b>Source</b>	Yeast
<b>Target Names</b>	Rdx
<b>Protein Names</b>	Recommended name: Radixin Alternative name(s): ESP10
<b>Expression Region</b>	1-583
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Radixin is a cytoskeletal protein that may be important in linking actin to the plasma membrane. It is highly similar in sequence to both ezrin and moesin. The radixin gene has been localized by fluorescence in situ hybridization to 11q23. A truncated version representing a pseudogene (RDXP2) was assigned to Xp21.3. Another pseudogene that seemed to lack introns (RDXP1) was mapped to 11p by Southern and PCR analyses.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.