




Recombinant Danio rerio Hydroxylysine kinase (hykk)

Product Code	CSB-EP428180DIL
Relevance	Catalyzes the GTP-dependent phosphorylation of 5-hydroxy-L-lysine.
Abbreviation	Recombinant Zebrafish hykk protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8WFT6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSTKESKPNSHSHSQVTDVVKRLYGLTASVVRPLPSYDDQNFYVAPSEGGEFV LKVMNSADSENVAVIELQTQSMNFLHQRLPAQTALPTLTGQLMSLEEFDCGF GTQIYLVRLTYLPGTTIAKITCSPQILYDVGKMAATLDTVLLQMEHPNTRVLQR ERFIWKLTSIPLLNQYVHVMDGDPVQKIVKGVIEKYQVQVMPKLPLFRECINHG DFNDHNLVLPDGPSPKYKISGILDFADMSCGYFIFELAITIMYMMIESPTPLDVG GPVVAGWESVFPLNEAERDSLWLVLMCRFCQSLVLARYTVIQQPENEEYLMIT SRKGLRHLRSLWELGKDEVERRWFQSAQQFRQI
Research Area	Cell Biology
Source	E.coli
Target Names	hykk
Protein Names	agphd1
Expression Region	1-355aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Mol. Weight	60.5 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.