



Recombinant Human Homeobox protein Hox-D3 (HOXD3)

Product Code	CSB-MP010686HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31249
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLFEQGQQAL ELPECTMQKA AYYENPGLFG GYGYSKTTDT YGYSTPHQPYP PPPAAASSLD TDYPGSACSI QSSAPLRAPA HKGAELNGSC MRPGTGNSQG GGGGSQPPGL NSEQQPPQPP PPPPTLPPSS PTNPGGGVPA KKPKGGPNAS SSSATISKQI FPWMKESRQN SKQKNSCATA GESCEDKSPP GPASKRVRTA YTSAQLVELE KEFHFNRYLC RPRRVEMANL LNLTERQIKI WFQNRMMKYK KDQKAKGILH SPASQSPERS PPLGGAAGHV AYSGQLPPVP GLAYDAPSPP AFAKSQPNMY GLAAYTAPLS SCLPQQKRYA APEFEPHPMA SNGGGFASAN LQGSPVYVGG NFVESMAPAS GPVFNLGHLS HPSSASVDYS CAAQIPGNHH HGPCDPHPTY TDLSAHHSSQ GRLPEAPKLT HL
Source	Mammalian cell
Target Names	HOXD3
Protein Names	Recommended name: Homeobox protein Hox-D3 Alternative name(s): Homeobox protein Hox-4A
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, located on different chromosomes, consisting of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXD genes located at 2q31-2q37 chromosome regions. Deletions that removed the entire HOXD gene cluster or 5 end of this cluster have been associated with severe limb and genital abnormalities. This protein may play a role in the regulation of cell adhesion processes.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.