




Recombinant Human 40S ribosomal protein S11 (RPS11)

Product Code	CSB-EP020365HU
Abbreviation	Recombinant Human RPS11 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P62280
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ADIQTERAYQKQPTIFQNKRVLLGETGKEKLPRYYKNIGLGFKTPKEAIEGTYI DKKCPFTGNVSIRGRILSGVVTMKMQRTIVIRRDYLHYIRKYNRFEKRHKNMMS VHLSPCFRDVQIGDIVTVGECRPLSKTVRFNVLVTKAAGTKKQFQKF
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	RPS11
Expression Region	2-158aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	45.3kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.