



Recombinant human Vesicular, overexpressed in cancer, prosurvival protein 1 protein

Product Code	CSB-MP035444h
Relevance	Increases the transcriptional activity of NFKB1 by facilitating its nuclear translocation, DNA-binding and associated apoptotic response, when overexpressed.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96AW1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	KKHCWYFEGLYPTYIICRSYEDCCGSRCCVRALSIQRLWYFWFLLMMGVLFCCGAGFFIRRRMYPPPLIEEPAFNVSYTRQPPNPGPGAQQPGPPYYTDPGGPGMNPVGNMAMAFQVPPNSPQGSVACPPPPAYCNTPPPPYEQVVKAK
Research Area	Epigenetics and Nuclear Signaling
Source	Mammalian cell
Target Names	VOPP1
Expression Region	23-172aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.